ES Manual de instrucciones
EN Instruction manual
FR Mode d'emploi
DE Bedienungsanleitung
IT Manuale d'istruzioni

WKE005

MELTY

CALENTADOR DE PARAFINA
PARAFFIN WARMER
CHAUFFE-PARAFFINE
PARAFFINBAD
SCALDAPARAFFINA

USO DEL MANUAL

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USE OF THE MANUAL

It is important that you read and retain this user manual for future reference. It contains important information on its use and maintenance.

You can also ask a professional for guidance to make sure you are using the bed correctly.

WARNING

- Never under any circumstances attempt to open or inspect the internal components of the device. If an inspection becomes necessary, please contact your supplier or a suitably competent technician.
- Never use the device in humid or poorly ventilated environments, such as wet rooms, bathrooms etc. to prevent damage to the electrical components.
- Never use the unit outdoors; exposure to the elements can cause damage.
- If the power cord was destroyed, please ask a professional to service and replace it to avoid hazards or contact the dealer for service.

EC-DECLARATION OF COMPLIANCE

This product meets the following directives:

- Low Voltage Directive: 2014/35/EU
- Electromagnetic Compatibility Directive: 2014/30/EU

RECYCLING

Whenever any of our electrical products or their presentation packaging bear this symbol, it means that these products should not be treated as conventional household waste in Europe. To ensure proper treatment of this waste, dispose of it in accordance with local laws or as required for the disposal of electrical equipment. This initiative will help preserve natural resources and improve standards of environmental protection regarding the treatment of electrical waste.

WELKOME!

Thank you very much for your purchase. Please read carefully the instructions before using the device. We accept no responsibility for damages caused by misuse.

The company reserves the right to amend the details of the devices with no prior notice. If any mistakes are found in these instructions, please feel free to send us an e-mail informing us so that we can correct them.

CONTENT

On arrival of the package, check the presence of the following items:

- Luxury paraffin bath

GENERAL INFORMATION

Beauty paraffin treatment is an application used both in beauty salons and in household beauty. The ingredients of the paraffin wax are vitamins, collagen, anti-wrinkle and whitening extracts that are beneficial for both foot and hand. Paraffin wax treatment is also called: “crystal mask”.

Paraffin wax is a petroleum by-product that has excellent heat retention properties. The paraffin coating causes an increase in blood circulation which helps the pore to open. As a result, this facilitates perspiration acting as a mini-sauna which makes the moisturizers be absorbed more effectively by the skin. The dry and rough skin results in a soft, fresh and healthy appearance.

REMARKS AND WARNING

- o Put the oval-shaped plastic plate on the inner bottom of the basin to prevent the foot from overheating.
 - o Before applying the paraffin on the body, you should check if the temperature is proper on your hand.
 - o The paraffin should not be used on broken skin.
 - o Do not pour the remaining paraffin from the container after the process. You should switch off the equipment first and wait for the paraffin to get cold. Then turn on the equipment again and wait until the paraffin at the bottom starts melting. Then you will be able to remove the whole block of paraffin easily.
-

OPERATION

1. Put the paraffin into the warmer.
2. Turn the switch "hi" to melt the paraffin.
3. When the paraffin melts, turn on the switch "lo". Once the temperature decreases, check the temperature with your hand to see if it is proper to start the treatment. Make sure it is not as high as to become uncomfortable for the patient.

HAND CARE RECOMMENDATIONS

1. Wash your hands and remove all the accessories you wear.
 2. Use the finger to check the paraffin temperature.
 3. Extend the whole palm into the paraffin for 4-6 times shortly for a 2-3 seconds contact. Then you will get a thick film on your hand.
 4. Put on plastic gloves to make sure that you keep the temperature of your hands so as to help the skin to absorb the nutrients. Take the gloves off after 15-20 minutes.
 5. Slowly remove the paraffin film from the hands and clean them carefully. The skin will become smooth and delicate.
-

FOOT CARE RECOMMENDATIONS

1. Wash the feet previously to the treatment.
 2. Use the finger to check the paraffin temperature.
 3. Put the feet into the paraffin 4-6 times shortly for a 2-3 seconds contact. If the feet are too big to put them into the container, you should use the brush to apply the wax on the feet. Once the paraffin forms a film on the feet, use a plastic bag to cover the whole feet in order to keep the temperature for around 15 minutes.
 4. Slowly remove the paraffin film from the feet and clean them carefully. The foot skin will become smooth and delicate
-

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